

Lansinoh®



SMARTPUMP™

DOUBLE ELECTRIC BREAST PUMP

App Pairing Instructions



First you'll need to install the Lansinoh® Baby App, which is a free software application for your phone that connects or “pairs” your pump with your smartphone. Once installed, the Lansinoh® Baby App will guide you through the registration and pump connection process and connecting your pump with the app.

To Install the App on an iPhone®

In order to minimize data usage, it is recommended to have a wifi connection while downloading the app. Download the app directly to your iPhone® from the App Store. Search for Lansinoh® Baby to find the app.

OR

Using your iPhone's browser, go to www.Lansinoh.com/App and click on the **Download on the Apple Store** icon.

Click the **GET** button and then the **INSTALL** button to download to your phone.

Once the app has finished installing, click the **OPEN** button.

To Install the App on an Android™ Phone:

In order to minimize data usage, it is recommended to have a wifi connection while downloading the app. Go to Google Play™ and search for Lansinoh® Baby.

OR

Using your phone's browser, go to www.Lansinoh.com/App and click on the **Google Play™** icon.

Select Install and follow the on-screen instructions.

Registering Your App

The app will automatically take you to the Login page. If you are a first time user, press **Click here to register**, which will open up the first time registration page.

Complete the registration page by entering your first and last name, email address, and create a password.


You may also choose to enter your baby's name, birthdate, & gender at this point, however this information is not required. It can be added at a later time by selecting **My Profile** in the drop-down menu on the upper-left corner of your screen.

If you intend to track more than one child, click on the **Add another baby** button and enter information for your second child.

Once the registration page is complete, click **Save** and you will automatically be taken to the **Add & Register Your Pump page**.

Tip: The app will need access to your camera in order to pair your pump. Be sure to allow access when your phone asks.

Tip: In order to receive tips, encouragement and important update information, allow the app to send you notifications if asked.



First time registration page

If you are using a Lansinoh® Smartpump™ Double Electric Breastpump, you can connect your pump to the app to automatically track your pumping sessions. If you are not using the Lansinoh® Smartpump™ select **'Not now' – return to homescreen** to use the app manually. To connect your Lansinoh® Smartpump™, complete the following steps:

Connecting Your Smartpump™

You will need to allow the app to have access to your camera in order to pair your pump. Turn the base of the pump over and locate the QR code on the bottom of the pump.

The app will automatically activate your phone's camera in order to allow you to scan the QR code and connect your pump. Once the camera function appears, place the pump base upside down on a hard surface and align the QR code inside the box.

Once the code is aligned, the pump will automatically begin to connect to the app and you will see a notification that your pump is connected.



Aligning and scanning the QR code

Adding a Pumping Session

Once paired, you will automatically be taken to the home screen of the Lansinoh® Baby App to begin use.

Once you have paired your phone and your pump, to log a pumping session, ensure you have a power source, and simply press the power button on your pump. Your pump will automatically begin recording your pumping session. Your pump will stop tracking when you press the power button again to turn the pump off. Your pumping session will be recorded immediately if your app and pump are paired.

When the pump is turned on, a white Bluetooth® icon in the upper left corner of your screen will confirm a Bluetooth® connection.



Connected



Not Connected

After your pumping session is complete, manually enter the volume of milk pumped by using the sliders on the milk collection bottles shown on the pumping input screen.

Click the **Save** button to complete your session.



Helpful Hints & Troubleshooting

The next time you turn your pump on, your Lansinoh® Smartpump™ should automatically connect to the app.

If you lose the connection to your pump, you can reconnect at any time by selecting **My Profile** in the drop down menu in the upper left corner of your screen and select **Add Smartpump**.

To sync your pump to your app, it is best to have the pump plugged into an electrical outlet. Simply turn the pump off when your pumping session is complete. If your pump does not have batteries, leave the pump plugged in for 1 minute after turning it off to ensure the pump has enough power to complete the data transfer.

Your pump will not record pumping sessions that are under 2 minutes in duration.

In order to maintain connectivity between your pump and your phone, they must be in a range of 10 feet of each other.

If you have not synced your app and your pump in more than 5 days, you may need to sync your phone prior to recording a new pumping session in order to maintain accurate time and date recordings.

If you keep batteries in your pump, time and date recordings should maintain accuracy regardless of the length of time between pumping sessions.

If you are unable to allow the app to access your phone's camera, contact Lansinoh's customer service team at customerservice@Lansinoh.com.

For more helpful tips and trouble shooting, visit www.Lansinoh.com/App to learn more.

©2016 Lansinoh Laboratories, Inc. All Rights Reserved.

iPhone and iTunes Store are registered trademarks of Apple Inc.

Android and Google Play are trademarks of Google Inc.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Lansinoh Laboratories is under license. Other trademarks and trade names are those of their respective owners.